

Primary Cord Naïve CD4+ CD45RA+ Cells

Prepared by: David Randolph (drdrdr@uab.edu) at University of Alabama, Birmingham

Goal: To obtain large numbers of highly pure primary CD4+ CD45RO- CD25- cells from cord blood, and prepare them as naïve, Th0, and Th1 cells. This protocol combines two purification strategies, a negative selection step with Rosette-Sep Reagents followed by an additional negative selection step with magnetic goat anti-mouse beads, to give desired purity but with higher yields than can be obtained with FACS sorting and cells are “untouched” with no antibodies bound to their surface.

Reagents:

- 1) Rosette-sep CD4 reagent
- 2) Lymphocyte separation media
- 3) Mouse anti-human CD45RO, mouse anti-human CD25, and mouse anti-human CD161 antibodies (eBioscience or BD Bioscience)
- 4) Goat anti-mouse magnetic particles (Chemicell)
- 5) Miltenyi LS magnetic separation columns
- 6) PBE buffer (PBS + 0.5%BSA + 2 mM EDTA)
- 7) T Cell Media = RPMI 1640 + 10%FCS + Pen/Strep + 2mM L-glut + 10mM Hepes + 50 μ M 2- β ME
- 8) 2X Freezing media = 8ml FCS + 2ml DMSO
- 9) Human recombinant IL-12 (R&D Systems)

Step 1: Negative Selection of CD4 Cells with Rosette Sep

- 1) Obtain sample of cord blood anti-coagulated with 10u/ml heparin
- 2) Dilute 1:2 with PBE
- 3) Underlay with lymphocyte separation media
- 4) Spin at 1200g x 30 minutes at RT. Collect CD4 T cells at interface.
- 5) Wash in PBE and count.

Step 2: Further enrichment with Gt anti-mouse beads

- 1) Add 0.05ug of anti-CD25, anti-CD161, and anti-CD45RO per million cells
- 2) Incubate on ice 30 minutes
- 3) Wash once with PBE
- 4) Add 3ul per million of goat anti-mouse magnetic particles.
- 5) Incubate on ice 30 minutes.
- 6) Wash once with PBE.
- 7) Resuspend in PBE, then run through Miltenyi LS column on magnet.
- 8) Wash column 3 times with PBE, collecting the flow-through.
- 9) Count cells, spin down, and resuspend in T Cell Media.

Step 3: Stimulation for DNase I experiments

- 1) For Naïve Cells, take 10 million in T Cell Media. Add 2X Freezing media and freeze immediately.

- 2) For Naïve Cells with Short Stim, Take 10 million and stimulate in T cell media at 37° and 7% CO₂ with plate-bound anti-CD3 + anti-CD28 for 4 hours. Then freeze cells.
- 3) For Th0 Cells, take 10 million and culture with CD3/CD28 beads with no cytokines added (neutral conditions) in T cell media at 37° and 7% CO₂ for 3 days. Remove beads, collect cells, and freeze.
- 4) For Th1 Cells, take 10 million and culture with CD3/CD28 beads + 1 ng/ml IL-12 (Th1 conditions) in T cell media at 37° degrees and 7% CO₂ for 3 days. Remove beads, collect cells, and freeze.