

## **Human peripheral blood-derived erythroblasts (PBDE)**

### **1. Cell expansion and Differentiation**

20mL blood was collected and the buffy coat isolated using ficoll-Paque.

Cells were expanded for 2 weeks as follows:

Week 1: Use StemSpan as basal medium and add:

- 10<sup>-6</sup> M Hydrocortisone
- 300ng/3mL SCF
- 120ng/3mL Flt-3L
- 40ng/3mL IL-3
- 40ng/3mL BMP-4
- 8U/3mL EPO

Week 2: Use StemSpan as basal medium and add:

- 10<sup>-6</sup> M Hydrocortisone
- 120ng/3mL SCF
- 120ng/3mL IGF1
- 40ng/3mL IL-3
- 40ng/3mL BMP-4
- 10U/3mL EPO

Feed the cells every other day. (Try to keep the density at  $\sim 2-3 \times 10^5$  /mL)

After 14 days, collect the cells, wash twice with PBS, and perform the crosslinking reaction.

### **2. Crosslinking reaction.**

The cells are collected by centrifugation at 2000-2500 x g for 10 min and crosslinked in PBS containing 1% formaldehyde for 10 min at room temperature on oscillating platform shaker.

The reaction is stopped by adding glycine to a final concentration of 0.125 M.

After 5 min, cells are washed twice with ice-cold PBS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free) and centrifuged at 1,500-2,000 x g for 5 min at 4°C.

Supernatant is discarded, and cell pellets are flash frozen in liquid nitrogen and stored at -80°C.

Note: StemSpan can be purchased at <http://www.stemcell.com>.